



Electro Surgical Instruments

RB Medikal Instruments



Laryngoscope, Otoscope, Ophthalmoscope

PRODUCTS CATALOG

Diagnostic Instruments



www.rbmedikal.com

RBmedikal

Electro Surgical Instruments

RB Medikal Instruments was established with a mission to provide our valued clients with impeccable service in the field of advance and update, finest quality hand crafted Surgical Instruments, Dental Instruments. Align with our team of professionals; we always strive for a Full satisfaction of a customer which is necessary for long terms business relations. We are also continuously developing new products, upgrading our instruments in all possible ways through technological enhancements and investigating under improved methods.

We are a company that sets high standards. Our growth has come from offering an unparalleled variety of innovative, high-quality products and services resulting from the dedication of our employees.

Full satisfaction of a customer is necessary for long terms business relations. We are committed to establish a life long partnership with all of our customers. Our commitment to quality, based on our continuous monitoring at every stage of the manufacturing process has placed us at par with international standards, while maintaining very competitive prices.

Certifications





LED/ Halogen Bulb

MACINTOSH

Integrated flexible Laryngoscope
(Fiber Optic) Cold light Guide.

MACINTOSH	Fig	Working Length
RB - 1001-01	1	72 mm
RB - 1001-02	2	88 mm
RB - 1001-03	3	112 mm
RB - 1001-04	4	134 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth edges. Rounded tip and high quality surface finish large fiber blades offers more light no black mark fully autoclavable and easy to clean stain finish.

MACINTOSH

Flexible interchangeable blade
(fiber optic) cold light guide.

MACINTOSH	Fig	Working Length
RB - 1003-02	2	80 mm
RB - 1003-03	3	92 mm
RB - 1003-04	4	116 mm
RB - 1003-05	5	135 mm

FEATURES:

303/304 Surgical stainless steel
Smooth edges. Rounded tip and high quality surface finish removable light guides design allow easier cleaning and extends blade life provides with 4mm fiber optic bundles for light pass gives true tissue color.



LED/ Halogen Bulb

MACINTOSH

Flexible fix tube fiber optic cold light guide.

MACINTOSH	Fig	Working Length
RB - 1005-01	1	70 mm
RB - 1005-02	2	88 mm
RB - 1005-03	3	112 mm
RB - 1005-04	4	134 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and high quality surface finish provides with 3mm Fiber Optic
Offers no black marks fully autoclavable and easier to clean stain finish.



LED/ Halogen Bulb



LED/ Halogen Bulb

MACINTOSH

Flexible fix tube fiber optic cold light guide.

MACINTOSH	Fig	Working Length
- 1007-02	2	92 mm
- 1007-03	3	113 mm
- 1007-04	4	135 mm
- 1007-05	5	155 mm


FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and high quality surface finish. Three spring ball bearings gives secure fit and firm support supplied with LED or halogen bulb satin finish.



Integrated Laryngoscope (Fiber Optic) Cold Light Guide

MACINTOSH	Fig	Working Length
RB - 1009-01	1	72 mm
RB - 1009-02	2	90 mm
RB - 1009-03	3	112 mm
RB - 1009-04	4	135 mm



FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and high quality surface finish large Fiber bundles offers more light no black marks fully Autoclavable and easier to clean satin finish.



LED/ Halogen Bulb




LED/ Halogen Bulb

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. high quality surface finish provides with 3mm Fiber Optic offers no Black Maks fully Autoclavable and easier to clean satin finish.

Miller Fix Tube HBT (Fiber Optic) Cold Light Guide

MILLER	Fig	Working Length
RB - 1011-00	00	45 mm
RB - 1011-0	0	55 mm
RB - 1011-01	1	80 mm
RB - 1011-02	2	130 mm
RB - 1011-03	3	170 mm
RB - 1011-04	4	180 mm





LED/ Halogen Bulb

Macintosh Interchangeable (Fiber Optic) Cold Light Guide

MACINTOSH	Fig	Working Length
RB - 1013-01	1	80 mm
RB - 1013-02	2	92 mm
RB - 1013-03	3	116 mm
RB - 1013-04	4	135 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and high quality surface finish. Removable Light guides design allow easier cleaning and extends blade life provides with 4mm Fiber Optic bundles for light pass gives true tissue color satin finish.



Miller Fix tube HBT (Fiber Optic) Cold Light Guide.

MILLER	Fig	Working Length
- 1015-00	00	45 mm
- 1015-0	0	55 mm
- 1015-01	1	80 mm
- 1015-02	2	130 mm
- 1015-03	3	170 mm
- 1015-04	4	180 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. High quality surface finish. Provides with 4mm Fiber Optic Bundles for light pass gives true tissue color fully autoclavable and easier to clean satin finish



LED/ Halogen Bulb



Macintosh Fix Tube (Fiber Optic) Cold Light Guide

MACINTOSH	Fig	Working Length
RB - 1017-01	1	70 mm
RB - 1017-02	2	88 mm
RB - 1017-03	3	112 mm
RB - 1017-04	4	134 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and High quality surface finish. Provides with 3mm Fiber Optic Bundles for light pass gives true tissue color fully autoclavable and easier to clean satin finish.



LED/ Halogen Bulb



Miller Fix Tube (Fiber Optic) Cold Light Guide

MILLER	Fig	Working Length
RB - 1019-00	00	45 mm
RB - 1019-0	0	57 mm
RB - 1019-01	1	80 mm
RB - 1019-02	2	130 mm
RB - 1019-03	3	170 mm
RB - 1019-04	4	180 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. High quality surface finish. Provides with 3mm Fiber Optic Bundles for light pass gives true tissue color fully autoclavable and easier to clean satin finish.



LED/ Halogen Bulb





LED/ Halogen Bulb

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and high quality surface finish. Three Spring ball bearing gives secure fit and firm support supplied with LED or Halogen Bulb
Satin finish.

Macintosh Conventional

MACINTOSH	Fig	Working Length
RB- 1021-00	00	60 mm
RB- 1021-0	0	68 mm
RB- 1021-01	1	92 mm
RB- 1021-02	2	113 mm
RB- 1021-03	3	135 mm
RB- 1021-04	4	155 mm



Miller Conventional


MILLER	Fig	Working Length
RB - 1023-00	00	60 mm
RB - 1023-0	0	70 mm
RB - 1023-01	1	92 mm
RB - 1023-02	2	113 mm
RB - 1023-03	3	135 mm
RB - 1023-04	4	155 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. High quality surface finish.
Three spring ball bearing gives secure fit and firm support supplied with LED or Halogen Bulb
Satin finish.



LED/ Halogen Bulb




Macintosh Conventional Eco

MACINTOSH	Fig	Working Length
RB- 1025-01	1	80 mm
RB- 1025-02	2	92 mm
RB- 1025-03	3	116 mm
RB- 1025-04	4	135 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and High quality surface finish. Three spring ball bearings gives secure fit and firm support supplied with LED 2.5v 5000 Lux satin finish




Macintosh Conventional Eco

MILLER	Fig	Working Length
RB- 1027-00	00	45 mm
RB- 1027-0	0	57 mm
RB- 1027-01	1	80 mm
RB- 1027-02	2	130 mm
RB- 1027-03	3	170 mm
RB- 1027-04	4	180 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. High quality surface finish.
Three spring ball bearing gives secure fit and firm support supplied with LED 2.5v 5000 Lux satin finish






Macintosh F/O Fix Tube (GS 100)
Cold Light Guide

MACINTOSH	Fig	Working Length
RB - 1029-01	1	80 mm
RB - 1029-02	2	92 mm
RB - 1029-03	3	116 mm
RB - 1029-04	4	135 mm

FEATURES

303/304 Surgical Stainless Steel
Smooth Edges. Rounded tip and high quality surface finish. Provides with 3mm Fiber Optic offer no black marks fully autoclavable and easier to clean finish



LED/ Halogen Bulb

Single USE Laryngoscope

Macintosh

Flexible conventional blade

MACINTOSH	Fig	Working Length
RB - 1031-01	1	76 mm
RB - 1031-02	2	88 mm
RB - 1031-03	3	112 mm
RB - 1031-04	4	130 mm

FEATURES

Premium quality Surgical Stainless Steel
Spatula provides with unbreakable plastic base
Unique heal design prevents Blade from touching the handle clearly
Marked and packed single patient use
smooth edges, rounded tip and high quality surface finish low cost then reusable one.



LED/ Halogen Bulb



MACINTOSH

Flexible Blade Fiber Optic

MACINTOSH	Fig	Working Length
RB - 1033-01	1	76 mm
RB - 1033-02	2	88 mm
RB - 1033-03	3	112 mm
RB - 1033-04	4	130 mm

FEATURES

Premium quality surgical stainless steel Spatula provides with unbreakable plastic base unique heal design prevents blade form touching the handle clearly marked and packed single patient use Fiber Optic tube black jacket end lit fiber smooth edges, rounded tip and high quality surface finish low cost reusable one.

Flexible Conventional Blade

MACINTOSH	Fig	Working Length
RB - 1035-01	1	76 mm
RB - 1035-02	2	88 mm
RB - 1035-03	3	112 mm
RB - 1035-04	4	130 mm

FEATURES

Premium quality surgical stainless steel spatula provides with unbreakable plastic base unique heal design prevents blade from touching the handle clearly marked and packed single patient use smooth edges, rounded tip and high quality surface finish low cost reusable one



LED Bulb



Macintosh F/O Fix Single Use

MACINTOSH	Fig	Working Length
RB - 1037-01	1	82 mm
RB - 1037-02	2	92 mm
RB - 1037-03	3	116 mm
RB - 1037-04	4	135 mm

Features

Tip to heal 303/304 Surgical Stainless steel
Smooth edges, Rounded tip and high quality surface finish.

Three spring ball bearings gives secure fit and firm support Fiber Optic tube back jacket and Lit fiber, satin finish.



LED/Halogen Bulb



Macintosh Eco Single Use

MACINTOSH	Fig	Working Length
RB - 1039-01	1	82 mm
RB - 1039-02	2	92 mm
RB - 1039-03	3	116 mm
RB - 1039-04	4	135 mm

Features

Tip to heal 303/304 Surgical stainless steel
smooth edges, rounded tip and high quality surface finish.

Three spring ball bearings gives secure fit and firm support plastic sleeve 2.5v LED 5000 Lux
Brighness satin finish.





Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
RB- 1041-01	1	79 mm
RB- 1041-02	2	88 mm
RB- 1041-03	3	112 mm
RB- 1041-04	4	130 mm

FEATURES

Tip to Spatula 303/304 Surgical Stainless Steel
Heel silver smooth edges, Rounded tip and High quality surface finish.
Three spring ball bearings gives secure fit and firm support fiber optic tube black jacket and LIT Fiber, satin finish



Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
RB - 1043-01	1	76 mm
RB - 1043-02	2	88 mm
RB - 1043-03	3	112 mm
RB - 1043-04	4	130 mm

Features

Tip to Spatula 303/304 Surgical Stainless Steel
Heel silver smooth edges, Rounded tip and High quality surface finish.
Three spring ball bearings gives secure fit and firm support fiber optic Acrylic Rod satin finish






Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
RB - 1045-01	1	76 mm
RB - 1045-02	2	88 mm
RB - 1045-03	3	112 mm
RB - 1045-04	4	130 mm

Features

Tip to Spatula 303/304 Surgical Stainless steel
heel plastic (ABS) smooth edges, rounded tip and
high quality surface finish.

Two spring ball bearings gives secure fit and
firm support fiber Optic tube black jacket and
LIT Fiber, satin finish



Miller Single Use (Fiber Optic)


MILLER	Fig	Working Length
RB - 1047-00	00	45 mm
RB - 1047-0	0	55 mm
RB - 1047-01	1	80 mm
RB - 1047-02	2	130 mm
RB - 1047-03	3	170 mm
RB - 1047-04	4	180 mm

Features

Spatula 303/304 Surgical Stainless steel
Heel plastic (ABS) smooth edges,
High quality surface finish

Two spring ball bearings gives secure fit and
firm support fiber optic tube black jacket end
Lit fiber , satin finish






Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
RB - 1049-01	1	76 mm
RB - 1049-02	2	88 mm
RB - 1049-03	3	112 mm
RB - 1049-04	4	130 mm

FEATURES

Tip to Spatula 303/304 Surgical Stainless Steel
Heel plastic (ABS) smooth edges, rounded tip and high quality surface finish.

Two Spring ball bearings gives secure fit and firm support fiber acrylic rod satin finish



Macintosh Conventional Single Use


MACINTOSH	Fig	Working Length
RB - 1051-01	1	83 mm
RB - 1051-02	2	95 mm
RB - 1051-03	3	120 mm
RB - 1051-04	4	138 mm

Features

Tip to Spatula 303/304 Surgical Stainless Steel
Unique Heel design smooth edges, rounded tip and high quality surface finish.

Two spring ball bearings gives secure fit and firm support 2.5 single and double LED Design satin finish





Miller Conventional Single Use

MILLER	Fig	Working Length
RB - 1053-00	00	45 mm
RB - 1053-0	0	55 mm
RB - 1053-01	1	82 mm
RB - 1053-02	2	130 mm
RB - 1053-03	3	135 mm
RB - 1053-04	4	185 mm

Features

Tip to Spatula 303/304 Surgical Stainless Steel
 Unique heel design smooth edges,
 High Quality surface finish
 Two spring ball bearings gives secure fit and firm support 2.5 single and double LED Design satin finish



Macintosh Conventional Single Use

MACINTOSH	Fig	Working Length
RB - 1055-01	1	80 mm
RB - 1055-02	2	92 mm
RB - 1055-03	3	116 mm
RB - 1055-04	4	133 mm

Features

Tip to Spatula 303/304 Surgical Stainless Steel
 Unique Heel Design Smooth Edges, Rounded tip and
 High Quality Surface Finish
 Two Spring ball hearings gives secure fit and firm support 2.5 Single LED Design 5000 Lux Brightness
 Satin Finish



HANDLES



Conventional Standard



Medium

Small

Stubby

RB - 1057-01

Fiber Standard



Medium

Small

Stubby

RB - 1057-02

Single Use Handles



Conventional

Conventional

Fiber Optic

RB - 1057-03

Hard Box & Soft Pouches



RB - 1059-01




RB - 1059-02



RB - 1059-03



RB - 1059-04



RB - 1060-01



RB - 1061-01



RB - 1061-02



RB - 1062-01



RB - 1063-01



RB - 1063-02



FEATURES

Equipped with an easy to use Bayonet Connecting system with high quality Professional metal head and chrome plated knurled battery handle making the surface a secure grip in the hand with extra bright light of 2.5v LED bulb illumination provided through Fiber Optic system. Three times high quality magnifying lens.

Requires: 20xC Batteries (Not included with set)
Packing: In Black soft zipper pouch with 10 specula having Dia 2.5, & 4.0mm

Fiber Optic
Conventional Otoscope
RB-1065-01



FEATURES

Equipped with an easy to use Bayonet Connecting system with high quality Professional metal head and chrome plated knurled battery handle making the surface a secure grip in the hand with extra bright light of 2.5v LED bulb illumination provided through Fiber Optic system. Three times high quality magnifying lens.

Requires: 20xC Batteries (Not included with set)
Packing: In Black soft zipper pouch with 10 specula having Dia 2.5, & 4.0mm

Fiber Optic
Metal Otoscope
RB-1066-01



FEATURES

Standard Metal Otoscope with bayonet connecting head fitting ideal for pneumatic tympanic membrane test and very day ENT use. 2.5v LED white light. Available in satin and polish finish.

Requires: 2xC Batteries (Not Included with set)
Packing: in black soft zipper pouch with 10 specula Having Dia 2.5mm & 4.0mm

Conventional
Metal Otoscope
RB-1067-01



FEATURES

High Quality Professional Parker Otoscope with extra bright light 2.5v LED Illumination Three times high quality magnifying lens. Available in satin and polish finish.

Requires: 2xC Batteries (Not included with set)
Packing: in black soft zipper pouch with 10 specula having Dia 2.5mm & 4.0mm

Conventional
Parker Metal Otoscope
RB-1068-01



FEATURES

High Quality Professional Parker Otoscope with screw fitting illumination light.
20.5v LED Metal Handle & Plastic head CPP.
Rheostat dimmer
Requires: 2xC Batteries (Not included with set)
Packing: In black soft zipper pouches with 10 Specula having Dia 2.5mm & 4.0mm

Conventional
Parker Metal Otoscope
RB-1070-01



FEATURES

2.5v Standard illumination LED Bulb.
Chrome Plated Handle and 4x Magnifying Glasses. Reusable Speculum Size 0, 2.5, 3.5 and 4mm include. Connecting for pneumatic otoscope battery handle with rheostat for adjustment of illumination.
Supply: German bulb technology for super view with 3 reusable ear specula.
Requires: 2xC Batteries (Not included with set)
Packing: In black soft zipper pouches
Available in satin and polish finish.

Conventional
Metal Otoscope
RB-1071-01



FEATURES

Open System with rotatable lens and optional speculum for convenient procedures
2.5v Otoscope.
Requires: 2xC Batteries (Not included with set)
Packing: In black soft zipper pouches with 3 VET reusable specula having Dia 4.0, 5.0, 7.0mm

Conventional
Veterinary Metal Otoscope
RB-1072-01



FEATURES

Having special focal distance with two fold magnification. Bayonet closure. 2.5v Otoscope.
Requires: 2xC Batteries (Not included with set)
Packing: In black soft zipper pouches with 3 VET reusable specula having Dia 4.0, 5.0, 7.0mm

Conventional
Veterinary Metal Otoscope
RB-1073-01



FEATURES

High Quality Pocket Sized mini Fiber Optic Otoscope with extra bright light 2.5v with 3x 4mm bundle with enhanced LED technology.

Requires: 2x AA Batteries (Not included with set)

Packing: In Gray or black pouch with 10 disposable specula Dia2.5mm & 40.0mm

Available Colors:

**Fiber Optic
Pocket Style**
RB-1075-01



**Fiber Optic
Pocket Style**
RB-1076-01



FEATURES

High Quality Pocket Sized mini Fiber Optic Otoscope with extra bright light 2.5v with 3x Magnification lens. Distal fiber optic illumination 4mm bundle with enhanced LED technology.

Requires: 2x AA Batteries (Not included with set)

Packing: In Gray or black pouch with 10 disposable specula Dia2.5mm & 40.0mm

Available Colors:

**Conventional
Pocket Style**
RB-1077-01



**Conventional
Pocket Style**
RB-1078-01



Ophthalmoscope
RB-1080-01

FEATURES

Top Quality Ophthalmoscope with fixation star aperture in pocket format. Quality optic in mini-format, Ideal for rounds and house calls. Attachment clip with integrated on/off switch. Secure, switches off automatically when returned to the pocket. High quality handle. white upper section/refined plastic shockproof, sturdy, non-slip. Enhanced LED Light Technology concentrated white light for perfect illumination and a brilliant image. 4 different apertures, A full-functin instrument incl, Fixation star, Range of lense, 20D to +20D, 2-piece handle and head, Easy to maintain, flexible

Replaceable Batteries Size AA
Supply: Bristle and Hard Plastic case Packing



Ophthalmoscope
RB-1081-01

The Ophthalmoscope is a tool used in medicine to examine the interior of the eye, studying the inner parts of the eye, such as the optic disk, retinal blood cessel, retinal, choroid and macula, can help diagnose diseases of the eyes.



RB-1082-01
FEATURES

GM Ophthalmoscope set pin contacting fitting ophthalmoscope head pin contact fitting lens selection wheel permits choice of corrective lenses from 0 to-25 and 0 to +40 diameters single large circle aperture, rheostat for adjustment of illumination 2 battery handle choice

Require: 2xC Size batteries (Not Included with sets)
Supply: Hard Case



RB-1083-01
FEATURES

Parker mini ophthalmoscope 2.5v LED Bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses -25 to +40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching for eye-glassess illumination lens dial for clear identification for adopter settings.

Require: 2 x AA Size batteries (Not Included with sets)
Supply: Soft pouch or plastic hard case.



RB-1084-01
FEATURES



Parker mini ophthalmoscope 2.5v LED Bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses -25 to +40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching for eye-glassess illumination lens dial for clear identification for adopter settings.



Require: 2 x C Size batteries (Not Included with sets)
Supply: Soft pouch or plastic hard case.




RB Medikal Instruments

 S.I.E, Fateh Garh, Near Umer Town, Sialkot- 51310 Pakistan

 +92 52 3552898  +92 52 3 576055

  +92-333-8619610

 rbmedikal@brain.net.pk / rbmedikal@gmail.com

 www.rbmedikal.com



RBmedikal

Electro Surgical Instruments



Offering the prime and astonishing collection of Diagnostic Instruments

www.rbmedikal.com

